

Chart 1: PCR Strip Tubes & 96-Well Plates Compatibility

●: Compatible ★: Optimal
Vc: Needs to be cut **Vn**: Compatible, but not

REAL-TIME PCR THERMAL CYCLERS		Product	PCR 8-Tube Strip				96-Well Non-Skirted Plate				96-Well Half-Skirted Plate				
			Cat. No.	F602001		F602002		F603101		F603102		F603105		F603106	
				0.1 ml		0.2 ml		0.1 ml		0.2 ml		0.1 ml		0.2 ml	
				Clear	White	Clear	White	Clear	White	Clear	White	Clear	White	Clear	White
Abbot	M2000rt			●	●										
ABI	StepOnePlus™ qPCR system	●	●			●	●								
	StepOne™ qPCR system (48-well)	●	●			Vc	Vc								
	ViiA™7 fast qPCR 96 x 0.1 ml	●	●			●	●								
	ViiA™7 qPCR system			●	●			●	●						
	7500 Fast qPCR system	●	●			●	●								
	7900HT Fast qPCR system	●	●			●	●								
	7500 qPCR system			●	●			●	●			●	●		
	7900HT qPCR system														
	QuantStudio® Q3/Q5/Q6/Q7			●	●			●	●						
	QuantStudio™ 12k Flex Fast	●	●			●	●								
	QuantStudio™ 12k Flex			●	●			●	●						
7300 qPCR system			●	●			●	●			●	●			
Agilent/Stratagene	Mx3000P™			●	●			●	●						
	Mx3005P™			●	●			●	●						
	Mx4000™			●	●			●	●						
Analytik-Jena	Toptical		●					●							
	qTower		●					●							
Bioer	Line-Gene I	Vn		●				●							
	Line-Gene K 48-well	Vn		●				●							
	Line-Gene 3310/3320	Vn		●				●							
	Line-Gene 9600	Vn		●				●							
BioRad	CFX96 Real-Time PCR	●	●			●	★			●	★				
	iQ5 Real-Time PCR			●	●								●		
	iQ Real-Time PCR			●	●								●		
	MyiQ Real-Time PCR			●	●								●		
	Opticon 2	●	●			●	●			●	●				
	CFX connect	●	★	●	●	●	★	●	●	●	●	●			
Eppendorf	Mastercycler ep realplex			●	●			●	●						
Hongshi Medical	SLAN-48P (48-well)			●				●							
	SLAN-96P (2 x 48-well)			●				●							
	SLAN-96S			●				●							
Roche	Cobas z 480	Vn	●								●				
	LightCycler® 480 System 384-well	Vn	●								●				
	LightCycler® 480 System	Vn	●								●				
	LightCycler® 96 System	Vn	●								●				
	LightCycler® Nano	Vn	●								●				
Tianlong Science	Gentier 96E			●	●			●	●						

STANDARD PCR THERMAL CYCLERS		Product	PCR 8-Tube Strip				96-Well Non-Skirted Plate				96-Well Half-Skirted Plate				
			Cat. No.	F602001		F602002		F603101		F603102		F603105		F603106	
				0.1 ml		0.2 ml		0.1 ml		0.2 ml		0.1 ml		0.2 ml	
				Clear	White	Clear	White	Clear	White	Clear	White	Clear	White	Clear	White
ABI	7000 ABI PRISM 7000 SDS 9700			●	●					●	●			●	●

Chart 2: PCR Strip Tubes Compatibility

●: Compatible

		Product	PCR 8-Tube Strip
		Cat. No.	F603001
		Volume	0.2 ml
		Color	Clear
ABI	StepOnePlus™ qPCR system		
	StepOne™ qPCR system		
	ViiA™7 fast qPCR 96 x 0.1 ml		
	ViiA™7 qPCR system		●
	7500 Fast qPCR system		
	7900HT Fast qPCR system		
	7500 qPCR system		●
	7900HT qPCR system		
	QuantStudio™ 12k Flex Fast		
	QuantStudio™ 12k Flex		●
	7300 qPCR system		●
	7000 ABI PRISM 7000 SDS		●
	9700		
	QuantStudio® Q3/Q5/Q6/Q7		●
	QuantStudio® Q3/Q5/Q6/Q7 fast		
	Veriti 96-Well/Veriti™ Dx 96-Well		●
	Veriti FAST 96-Well/Veriti™ Dx FAST		
SimpliAmp™		●	
Agilent/Stratagene	Mx3000P™		●
	Mx3005P™		●
	Mx4000™		●
	Agilent AriaMX		
	SureCycler 8800 thermal cycler		●
Analytik-Jena	Toptical		
	qTower		●
	Tone96/96G		●
	Alpha SC		●
	FlexCycler		●
	Alpha SC		●
	TRIO 48		●
	TRIO Combi		●
SpeedCycler 36		●	
Bioer	Line-Gene I		●
	Line-Gene K 48-well		●
	Line-Gene 3310/3320		●
	Line-Gene 9600		●
	GenePro B		●
	Life ECO Thermal Cycler		●
	Life Express		●
	LifePro TC-96/G/H(b)A		●
	LifeTouch TC-96/G/H(b)B		●
XPcycler TC-XP		●	
BioRad	CFX96 Real-Time PCR		
	iQ5 Real-Time PCR		●
	iQ Real-Time PCR		●
	MyiQ Real-Time PCR		●
	MyiQ2 Real-Time PCR		●
	Opticon 2		
	CFX connect		●
	T100		●
	C1000/ C1000 Touch		●
	DNA Engine		●
	MiniOpticon		●
	MJ Mini		●
	PTC-100		●
	S1000		●
iCycler		●	
Eppendorf	Mastercycler ep realplex		●
	Mastercycler nexus /nexus eco /nexus GSX1		●
	Mastercycler ep gradient		●
	Mastercycler pro/pro S		●
	Thermomixer C		●
	Thermomixer comfort		●
Mastercycler gradient		●	
Hongshi Medical	SLAN-48P (48-well)		●
	SLAN-96P (2 x 48-well)		●
	SLAN-96S		●
Roche	Cobas z 480		
	LightCycler® 480 System 384-well		
	LightCycler® 480 System		
	LightCycler® 96 System		
	LightCycler® Nano		

Chart 3: PCR 384-Well Plates Compatibility

●: Compatible

		Product	384-Well Skirted Plate	
		Cat. No.	F600384	F601384
		Volume	40 µl	
		Color	White	Clear
ABI	7900HT/Fast			●
	StepOnePlus™			●
	ViiA™ 7/Fast			●
	3100			●
	3130			●
	3700			●
	3730			●
	3730x			●
	9700			●
	Veriti®/Fast			●
Agilent	SureCycler® 8800			●
Amersham/GE	MegaBACE 4000			●
Analytik-Jena/Biometra	FlexCycler			●
	TProfessional Standard/Basic/TRIO			●
	TRobot			●
	Uno II			●
BioRad/MJ Research	C1000 Touch™			●
	DNA Engine Dyad®/Dyad Disciple™			●
	Engine Tetrad® 2			●
	PTC-100®			●
	PTC-200			●
	PTC-225 Tetrad			●
	S1000™			●
CFX384™ Touch			●	
Corbett Life Science	Palm-Cycler™ 384			●
Eppendorf	Mastercycler M384			●
	Mastercycler pro 384			●
	Mastercycler pro/S			●
MWG	Primus 384			●
PeqLab	peqSTAR 384			●
Qiagen	FlexiGene®			●
Roche	LightCycler® 480 System 384-well	●		
SensoQuest	Thermocycler 384			●
	Thermocycler 96			●
Techne	Genius			●
	Genius Quad			●
	Prime			●
	Prime Elite			●
	Prime Elite Satellite			●
	PrimeG			●
	Touchgene® Gradient (TC512)			●
	Touchgene® X			●
Thermo Hybaid	Multiblock System			●
	Omnigene			●
	PCR Express			●
	Px2			●
	PxE			●
	Touchdown			●